



ONeal Medical - CPAP Plus

Regarding Minimum Pressure Support for Auto Bilevel

We often see auto bilevel prescriptions that request a ‘minimum pressure support.’ However, ONeal Medical’s supplier of auto bilevel machines, the Resmed AirCurve10 VAuto, *does not require minimum pressure support to be defined*; instead, the device automatically adjusts it based on Resmed’s proprietary algorithm. Thus, a prescribed minimum pressure support is unnecessary for the Resmed AirCurve10 VAuto.

With Resmed, maximum pressure support is set to patient comfort, with a default setting starting at 4 cm H₂O (or greater).

Below is an excerpt from page 10 of the Resmed Script Conversion Guide, which clarifies the difference between settings commonly prescribed for Philips and Resmed machines. The entirety of this document is available at onealmed.com/sleep.

Scripted therapy mode: Auto-bilevel

| Mode | Philips | ResMed recommended setting | | Mode |
|--------------|-------------------------------------|-------------------------------------|---|---|
| Auto-bilevel | DreamStation Auto BiPAP | AirCurve 10 VAuto | Conversion formula | VAuto <small>Note: CPAP and S modes also available on device</small> |
| | IPAP Max (4–25 cm H ₂ O) | Max IPAP (4–25 cm H ₂ O) | Max IPAP = IPAP Max | |
| | EPAP Min (4–25 cm H ₂ O) | Min EPAP (4–25 cm H ₂ O) | Min EPAP = EPAP Min | |
| | PS Min (0–8 cm H ₂ O) | | | |
| | PS Max (0–8 cm H ₂ O) | PS (0–10 cm H ₂ O) | Adjust PS for comfort; selected PS should not exceed Max PS and should ideally be 4 cm H ₂ O or greater to provide expiratory relief | |
| | Flex (Bi-Flex, None) | | | |
| | Flex Level (1, 2, 3) | | | |

Example

| | | | | |
|--------------|----------------------------------|----------------------------------|---|-------|
| Auto-bilevel | DreamStation Auto BiPAP | AirCurve 10 VAuto | Conversion formula | VAuto |
| | IPAP Max: 25 cm H ₂ O | Max IPAP: 25 cm H ₂ O | Max IPAP = IPAP Max | |
| | EPAP Min: 6 cm H ₂ O | Min EPAP: 6 cm H ₂ O | Min EPAP = EPAP Min | |
| | PS Min: 0 | | | |
| | PS Max: 6 | PS: 4 cm H ₂ O | Adjust PS for comfort; selected PS should not exceed Max PS and should ideally be 4 cm H ₂ O or greater to provide expiratory relief | |
| | Flex: Bi-Flex | | | |
| | Flex Level: 2 | | | |