

# Script note pad

**Patient name:** \_\_\_\_\_ **Order date:** \_\_\_\_\_  
**Address:** \_\_\_\_\_ **DOB:** \_\_\_\_\_  
**Email:** \_\_\_\_\_ **Phone number:** \_\_\_\_\_  
**Insurance name:** \_\_\_\_\_ **Insurance number:** \_\_\_\_\_  
**Date seen prior to sleep test:** \_\_\_\_\_ **Diagnosis test date:** \_\_\_\_\_  
**Duration of need:**  99 mo.  Other: \_\_\_\_\_

**Diagnosis:**  OSA: ICD-10 G47.33  Primary CSA: ICD-10 G47.31  
 Nonobstructive alveolar hypoventilation: ICD-10 G47.34  Other: \_\_\_\_\_

**CPAP therapy E0601 + E0562 + A9279**

**AirSense™ 10 Elite** with heated humidification

Pressure: \_\_\_\_\_ cm H<sub>2</sub>O (4–20 cm H<sub>2</sub>O) **DEFAULTS 10**  
 Ramp time: \_\_\_\_\_ min(s) (Auto, OFF–45 min.)  
 EPR™:  1  2  3

**APAP therapy E0601 + E0562 + A9279**

**AirSense™ 10 AutoSet™** with heated humidification

AutoSet mode **DEFAULTS**  
 Default mode settings  
 Min. pressure: \_\_\_\_\_ cm H<sub>2</sub>O (4–20 cm H<sub>2</sub>O) **4**  
 Max. pressure: \_\_\_\_\_ cm H<sub>2</sub>O (4–20 cm H<sub>2</sub>O) **20**  
 Ramp time: \_\_\_\_\_ min(s) (Auto, OFF–45 min.)  
 EPR™:  1  2  3

**AirSense™ 10 AutoSet™ for Her** with heated humidification

AutoSet for Her mode **DEFAULTS**  
 Default mode settings  
 Min. pressure: \_\_\_\_\_ cm H<sub>2</sub>O (4–20 cm H<sub>2</sub>O) **4**  
 Max. pressure: \_\_\_\_\_ cm H<sub>2</sub>O (4–20 cm H<sub>2</sub>O) **20**  
 Ramp time: \_\_\_\_\_ min(s) (Auto, OFF–45 min.)  
 EPR™:  1  2  3

**Bilevel therapy E0470 + E0562 + A9279**

**AirCurve™ 10 VAuto** with heated humidification

VAuto mode **DEFAULTS**  
 Default mode settings  
 Max. IPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–25 cm H<sub>2</sub>O) **25**  
 Min. EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–25 cm H<sub>2</sub>O) **4**  
 PS: \_\_\_\_\_ cm H<sub>2</sub>O (0–10 cm H<sub>2</sub>O) **4**  
 Ramp time: \_\_\_\_\_ min(s) (OFF–45 min.)

Spont mode  
 IPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–25 cm H<sub>2</sub>O) **10**  
 EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (3–25 cm H<sub>2</sub>O) **4**  
 Ramp time: \_\_\_\_\_ min(s) (OFF–45 min.)  
 Easy-Breathe ON

**AirCurve™ 10 S** with heated humidification

Spont mode  
 IPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–25 cm H<sub>2</sub>O) **10**  
 EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (3–25 cm H<sub>2</sub>O) **4**  
 Ramp time: \_\_\_\_\_ min(s) (OFF–45 min.)  
 Easy-Breathe ON

**S9™ VPAP™ COPD\*** with heated humidification

Spont mode **DEFAULTS**  
 Default mode settings  
 IPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–30 cm H<sub>2</sub>O) **13**  
 EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (3–25 cm H<sub>2</sub>O) **5**  
 Ramp time: \_\_\_\_\_ min(s) (OFF–45 min.)

**Bilevel w/ backup rate therapy E0471 + E0562 + A9279**

**AirCurve™ 10 ASV** with heated humidification

ASV mode **DEFAULTS**  
 Default mode settings  
 EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–15 cm H<sub>2</sub>O) **5**  
 Min. PS: \_\_\_\_\_ cm H<sub>2</sub>O (0–6 cm H<sub>2</sub>O) **3**  
 Max. PS: \_\_\_\_\_ cm H<sub>2</sub>O (5–20 cm H<sub>2</sub>O) **15**  
 Ramp time: \_\_\_\_\_ min(s) (OFF–45 min.)  
 Backup rate: Automatic (15 BPM)

ASV Auto mode **DEFAULTS**  
 Default mode settings  
 Min. EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–15 cm H<sub>2</sub>O) **4**  
 Max. EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–15 cm H<sub>2</sub>O) **15**  
 Min. PS: \_\_\_\_\_ cm H<sub>2</sub>O (0–6 cm H<sub>2</sub>O) **3**  
 Max. PS: \_\_\_\_\_ cm H<sub>2</sub>O (5–20 cm H<sub>2</sub>O) **15**  
 Ramp time: \_\_\_\_\_ min(s) (OFF–45 min.)  
 Backup rate: Automatic (15 BPM)

**AirCurve™ 10 ST** with heated humidification

Spont/timed mode  
 IPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–25 cm H<sub>2</sub>O) **10**  
 EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (3–25 cm H<sub>2</sub>O) **4**  
 Rate: \_\_\_\_\_ BPM (5–50 BPM) **10**

**AirCurve™ 10 ST-A** with heated humidification

Timed mode  PAC mode

Spont/timed mode  
 IPAP: \_\_\_\_\_ cm H<sub>2</sub>O (4–30 cm H<sub>2</sub>O) **10**  
 EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (3–25 cm H<sub>2</sub>O) **4**  
 Rate: \_\_\_\_\_ BPM (5–50 BPM) **10**  
 Ti: \_\_\_\_\_ sec. (0.1–4 sec.) **2**

iVAPS™ mode  
 Height: \_\_\_\_\_ in. (44–100 in.) **70**  
 Target patient rate: \_\_\_\_\_ BPM (8–30 BPM) **15**  
 Target Va: \_\_\_\_\_ L/min. (1–30 L/min.) **5.2**  
 Vt (Tidal volume) \_\_\_\_\_ (mL)  
 Vt/kg \_\_\_\_\_ (mL/kg)  
 EPAP: \_\_\_\_\_ cm H<sub>2</sub>O (3–25 cm H<sub>2</sub>O) **4**  
 Min. PS: \_\_\_\_\_ cm H<sub>2</sub>O (0–20 cm H<sub>2</sub>O) **4**  
 Max. PS: \_\_\_\_\_ cm H<sub>2</sub>O (0–27 cm H<sub>2</sub>O) **20**

\*Device does not include integrated wireless module, but does include SD card.